



Nasopharyngoscope

Excellent image meets flexibility

700715FX Nasopharyngoscope

- Insertion tube: Ø 3.4 mm
- Working length: 300 mm
- Field of view: 85°
- Angulation up/down: 160°/160°
- Standard accessories: Leakage tester, ETO-Cap

- **High-quality, German made**
- **Sharp image (18.000 pixels)**
- **Ergonomic design**



NEW for OLD! Get 25% discount on exchange scopes
(only valid for trade-in of western brand nasopharyngoscopes like Storz, Olympus,...)